



Instruction for use

Please read this information carefully because it will help you to use HYLO NIGHT® correctly

HYLO NIGHT®

For improvement of the tear film and protection of the eye's surface in cases of dry eye disease

Eye ointment without preservative

PRODUCT DESCRIPTION

HYLO NIGHT® is a sterile eye ointment with Vitamin A for ocular use.

COMPOSITION

1 g HYLO NIGHT® contains 250 I.U. Retinol Palmitate (Vitamin A), Liquid Paraffin, Light Liquid Paraffin, Wool Fat and White Soft Paraffin.

1 dosage = 1 strip of eye ointment (approx. 1 cm long) corresponds to approx. 5 I.U. Vitamin A.
One tube contains 5g of eye ointment.

INTENDED PURPOSE

HYLO NIGHT® improves the tear film and protects the eye's surface in cases of dry eye disease.

MEDICAL INDICATION

Dry Eye Disease

SUITABLE FOR

HYLO NIGHT® is suitable for patients of all ages. Based on existing knowledge, children, pregnant and breast-feeding women can also use HYLO NIGHT®.



BENEFITS, PERFORMANCE AND EFFECTIVENESS

The clinical benefit is the relief from symptoms of dry eye disease.

HYLO NIGHT® is a particularly soft and smooth eye ointment that spreads very well on the eye's surface.

Vitamin A is a natural component of the tear film supporting the eye ointment's good tolerance and helps the ointment to mix well with the available tears. Thus, the sensation of burning, dry or tired eyes is relieved. The eye lids slide easily over the cornea and conjunctiva.

Air conditioning, e.g. in cars or airplanes, wind, cold, direct sunlight, pollution or cigarette smoke often cause dry or itching eyes with a feeling of dryness or foreign bodies.

Working in front of a screen for long periods, watching TV or driving at night can also cause these symptoms. If you are wearing contact lenses, there is an even higher risk for these impairments of the eye's condition.

HYLO NIGHT® relieves these symptoms and constitutes a convenient long-lasting film on the eyes.

The fatty consistency of the ointment causes a slight impairment of vision. Therefore, HYLO NIGHT® should be applied before going to bed. Possible discomfort when opening the eyes for the first time in the morning can be prevented by using HYLO NIGHT®.

GENERAL INSTRUCTIONS FOR USE

Eye ointments should be used avoiding contact of the tube-tip with the eye or skin.

Only use HYLO NIGHT® if the seal of the carton was intact before the first use (☹).

Wash your hands and stay in front of a mirror.

Unscrew the cap covering the nozzle. Apply the eye ointment by gently squeezing the tube without folding or rolling.



Tilt your head backwards and gently pull on the lower eyelid using a finger of one hand.



With the other hand hold the tube in an upright position above the eye. Instill one strip of eye ointment into the gap formed between the eye and eyelid by gently squeezing on the tube.



Close your eye gently then blink normally.

Replace the cap firmly back onto the tube.

HYLO® NIGHT can be used in combination with lubricating eye drops, however it should not be used with other eye drops at the same time. HYLO NIGHT® should be used about 30 minutes after the application of the eye drops. If you need to use another eye ointment, please keep an interval of at least one hour between both applications.

HYLO NIGHT® can be used for 6 months after first opening and should not be used after the expiry date (📅). Please contact the manufacturer, your pharmacist or ophthalmologist if you have any questions about the correct use of this medical device.





DOSAGE

The frequency of use depends on your personal needs. HYLO NIGHT® should be used once a day before going to bed. In more severe dry eye cases, HYLO NIGHT® can be applied three times a day into each eye.

If you feel the need to use HYLO NIGHT® more often, speak your doctor.

Treatment can continue as long as needed.

If you follow the dosage instructions, overdosing with HYLO NIGHT® is unlikely as the target tissue is the ocular surface, especially the conjunctival sac, which has limited volume.

PRECAUTIONS AND WARNINGS

If your symptoms do not improve or if you feel worse after several days of use, talk to your doctor or pharmacist.

If you use HYLO NIGHT® during the day, your vision will be impaired for some time because of the consistency of the eye ointment. You should not drive a car, use machines or work in an unsecure environment during this time.

Only one person should be treated with one tube HYLO NIGHT® ().

SIDE EFFECTS

The ingredient Wool Fat may cause localised skin reactions (e.g. contact dermatitis). Please check before use if you have a known hypersensitivity to Wool Fat.

In very rare cases hypersensitivity reactions like redness, burning and/or itching eyes or excessive tearing and other temporary local irritations have been reported which improved after the treatment with HYLO NIGHT® was stopped.

If a serious incident has occurred during the use of this medical device or as a result of its use, please report it to the manufacturer () and to your national authority.

CONTRAINDICATIONS

HYLO NIGHT® should not be applied

- if you are hypersensitive to any of the ingredients
- while wearing contact lenses
- directly after ophthalmic-surgical interventions or in acute injuries of the eye as no clinical studies confirm this indication.

STORAGE

Do not store above 25°C ( 25 °C).

Keep HYLO NIGHT® out of the reach and sight of children. (imprinted on the packaging carton).

Any medical device you no longer use should be thrown away via household waste. This measure will help to protect the environment.

DATE OF INFORMATION:

July 2025

 URSAPHARM Arzneimittel GmbH
Industriestraße, D-66129 Saarbrücken
www.ursapharm.de

in UK by UK Responsible Person:

Scope Ophthalmics Ltd
1st Floor, Unit 4 Amberley Court
Whitworth Road, County Oak Way
Crawley, West Sussex, RH11 7XL
Freephone UK: 0808 175 3200
Fax: +44 203 515 9461

in Ireland

Scope Ophthalmics Ltd
Suite 5, Westland House,
Westland Park, Willow Road,
Dublin 12, Ireland
Telephone Ire: 1800 816 005
Email: info@scopeeyecare.com
Website: www.scopeeyecare.com

CE₀₁₉₇

MD

STERILE A

Legend of pictograms:

Symbol	Explanation
	Distributor
MD	Medical device
CE	CE mark
UDI	Unique Device Identifier printed on the folding carton.
STERILE A	Sterilized using aseptic processing techniques
LOT	Lot number
	Date of manufacture [YYYY/MM]
	Consult instructions for use

